



## Cell Cycle Staining Buffer

### I.

: CCS01 : 50ml/50T  
: 2 - 8°C :

### II.

(Propidium iodide, PI) DNA 4 - 5  
PI RNA DNA  
20 - 30 ( ) RNA DNA  
DNA PI PI stokes  
PI mRNA  
PI RNase A 50 ml DNA 50

### III.

1.  $2 \times 10^5 - 1 \times 10^6$  PBS
2. 1 ml DNA staining solution 5 - 10 30  
1

### IV.

- 1.
- 2.
- 3.

### V.

CCS01	Cell Cycle Staining Buffer	50 ml/50 T
CCS012	Cell Cycle Staining Kit	50T

	CCS01	CCS012
DNA staining solution (A)		
Permeabilization solution (B)		
		/

DNA CCS01 DNA CCS012  
CCS012 CCS012 CCS01

Product For Research Use Only rev: CCS01-01

MultiSciences Biotech Co., Ltd

Technical Service: tech@liankebio.com

To place an order: order@liankebio.com